

## LRRC37B-leucine rich repeat containing 37B Gene

Reference: PTXBC117373

|                |                                                                                                                          |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | LRRC37; leucine-rich repeat-containing protein 37B; C66 SLIT-like testicular protein; leucine rich repeat containing 37B |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | leucine rich repeat containing 37B                                                               |
| <b>GENE ID</b>          | 114659                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC117373                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | LRRC37B                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcttggctgctgttctggggcccatggccctctctacgtggcaactattgtcttactagtcaggaggtcagcctctgtgtgtgggtcaaggacccgtccagctgaccttaacccccctggggccacctgagccctgtgtctccctctccatctccatgggaatctccccatgc  
acctgtccccccagcagccccggggacttggattacctggggccctctgtctctcgcagatgtcagccctgctcagggaaccaactgaaaatttggctccattctgaagggaattggattcagctggagagctgccccctggggccagagccgttcttggctgcacatcaggacttaaatga  
caagcggactccagaagaaggtctccagaggtgtgtctcgtctcaaccgggacagaaacaggccctagtgcagctctcgcctcaagtggtgttcaactacagatctagatcgggctgcaggtcatcaggcagatgaaatcttgttccactagacagtaagtttcaagcaacca  
aatgtgtgttctgcccaagaacctgaaagaagatcagctgaaacgttggagcctctcgtgagattgttgggtccacaccaattatccaaacctcagcgtcagaagaacagacttggccagatgattatttggatgacacactgtatccggcagcctacctccagaactcgggtgaacgca  
gatgagcctccaggccctctgagcaagttggacttctcaattccatctagagcccaaaagtcacaaatccagagacccttgaagacatccagctctcttactccagggaagaagccccagcgcagcttctacagctccctcaggaggtagaaccttaacccagcaggagggcccgactc  
tgcctccagagctctctatggagagctcagctcaaacctccactgaatcatgaagtgacagttcaacctccaggtgaggtacagctcattataatttggcccaagtttaccagttcaaacctgcagatgtggaggttaccatgactcagagcctaaaaatgagacagaatctaccaagcccgca  
ggagggccccaattcagcctcccgaggagggcggaaccttctctacagccctgaggactacagatcctcctccagaacacctgagggtgacacttccaccttcagacaagggtcaggctcagcattcacacctgactgaagccacagttcaacctctggacctggagcttagcataactaca  
gagcctactacagaggttaaacctctccaaaccggagggaacctcagctcagcctccagaccgggggtgcccataactccagaaccactacagagattggacattccacacctggagaagactagagctctctacagaccaggttcagactctgcatcgaagcctgactgaa  
gtcacagggtccacctcaaaagttagaattcttcgaggattcattgtgtgagctgaaactgcaccagagggaacagaaagccctccacaagcaccacatatgtgagctctgacctcggagatgagactctgcatgtgtgtgtcagccccaaagcagaggtctcccaagtgctgtgc  
cagagcccgacacctacaatgcatcttcaccaccttaatttccaaggaaactatatttcataaccttgatgaaatgtatggaaagcatacagttggaccgagaaactaatttcagtgaaaattattgactgaattacctaaggttcaattgaaagcctgtataacctccagatttaattctcaa  
tcgcaatctctgactactgtcgaagatccatctcttgaactcggcgattaaatattctagacatgggaacaacacacatcacacttaacacacttaagaacattctcagatgactgttgaactggaaaaactgatcttaccagccatatggcctgtcctctcctccaatttaaaaatagca  
ttgaggtgtctgcaagacagtcagctgattgcaacactgcatgtctgactaacagcagataattgtctgaagaagcatctgtagggaatccagaaggagcgtttcagaagattgtacaagccggagcagcacatgagcactcagctgactattgagtcggagcgccctcagacag  
cagtggtcatcaactgtcaggtcttgggggtgatcagcttgaattcagctaacgcagctacgttccctcatcccaacgaggatgtgagaaggttcagctcaggttccggaccttgaaaatggaatgttcagaaacacatgtcgaagggagctgtgccaagctcatgttgcgaaca  
ggcctctgatgaagcttctcagcagcagcaggaagcaaaagcattgaatgtgaatgggatacggaccaacaaaaaacaattatattatgagaacatggaacagaatgaacagaaagcagaagtcagtgactgaagaaggttcaggagatgactataagaacaaact  
catcttcgcaatctctgactgtaataactaattttgattataattttgtcttataagggtgaattacataaaagggtcagagaataacaaagcaacccatcaatcaggagcctga

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)