

FLJ11292-hypothetical protein FLJ11292 Gene

Reference: PTXBC117434

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ11292
GENE ID	55338
ORIGIN SPECIES	Human
ACCESSIONS	BC117434
VECTOR	pDONR223
NCBI SYMBOL	FLJ11292
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggactctggcagcaatttcaaaactgtatgccagacctggcccttgcagtgatataaaattttcttgaggtgtctggaagatggctgaataggaaacagctctgtctgcagctccagcgagatcaacgcagaaggagataatttctgatttccaactaaagtaccagctcatctcatt
gggactgggttagacagtgggtgcagccacagaaggcaagcagaagcagggtgggtgtcgcctcaccgggaagcgaagggtcagggaactccctccctagccaagggaagctgtgaggactgtgaggaacggggcattccggcagacatactatgtttcccca
cgggtcttgcacccacagaccaggagattcccttgggtgctgcaccaccaggccctgggttcaagcaaaaaactgggcagccatttggcagacactga