

ZDHHC22-zinc finger, DHHC-type containing 22 Gene

Reference: PTXBC117676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	palmitoyltransferase ZDHHC22; C14orf59; DHHC-22; zinc finger DHHC domain-containing protein 22; zinc finger, DHHC domain containing 22; zinc finger DHHC-type containing 22		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 22
GENE ID	283576
ORIGIN SPECIES	Human
ACCESSIONS	BC117676
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtggccctgcggctgtcaacgtgtgtggccccgcctactttgtgcatctccctgtgtgaccttctgtgtgcagctcttctcttcttccctgccagcatgcgcgaggaccccgccgcccggctcttctgcgcgcctgtctccacggggcgctcttctattctctcggccaacgecc
tgggcaattacgtctgtcatccagaactccccagagacgtgggggctgtccagggggctcggccagggaagactccatgccctcactagcaccacttctgcccagtggtgtgcccagagtcacccctgaggcagaccatcactgtttcttcaccggcaactgcatcggcagcagg
aacatgcgcaactctgtctgttctgacctacacccctggcctgctctactccatgggtggccggcgtggcctacatctccgtgtcttcttccatctcttgcaccccttggccttctcagctctcgtccacatccatcagccagttcttctccggagctgtcctgggtctgaaatgttc
gtcatctcatgtctacctctgttgcgccatcgccctggcctgcgcggccttctgtgtccaccagctgtgtgtatctccggcggcagaccgcccaccaggtgcggaagggtgtgagtgagggcccgccctggcgaagaactacaagaggtcttcgaaagaggtggctgtct
gggctctgtgttcccatgttcaatgtcgaagtgagagctccaagcagcaggataagtag