



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LIN7A-lin-7 homolog A (C. elegans) Gene

Reference: PTXBC118609

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIN-7A; LIN7; MALS-1; TIP-33; VELI1; protein lin-7 homolog A; mammalian LIN-7 1; mammalian lin-seven protein 1; tax interaction protein 33; vertebrate LIN7 homolog 1; lin-7 homolog A, crumbs cell polarity complex component		

For research use only.

GENE DESCRIPTION	lin-7 homolog A (C. elegans)
GENE ID	8825
ORIGIN SPECIES	Human
ACCESSIONS	BC118609
VECTOR	pDONR223
NCBI SYMBOL	LIN7A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr