

LOC653166-olfactory receptor pseudogene Gene

Reference: PTXBC118668

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	olfactory receptor pseudogene
GENE ID	653166
ORIGIN SPECIES	Human
ACCESSIONS	BC118668
VECTOR	pDONR223
NCBI SYMBOL	LOC653166
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatggagataaccagagtgagaactcacagttcctctctgggatctcagagagtcagagcagcagatcctgttttgatgttcctgtccatgtacctgggtcacgggtctgggaatgtctcatcatcctggccatcagctctgattccacctgcacaccccatgtactcttc
ctggccaacctctccttcaactgacctctcttgcaccaacacaaatcccaagatgctgggaactccagtcaccagaacaaagccatctctatgcagggtgtctgacacagctctactctctgtctcttgggtgacctggacaacctcatctggcctgatggcgtatgatcgtatgtgg
ccatctgtgccccctccactatgtcacagccatgagccctgggctctgtgcttctctctctgtgttggggctgtctgtctctatggcctctctcactctctctgaccagggtgacctctgtgggccttga