

NCOA3-nuclear receptor coactivator 3 Gene

Reference: PTXBC119001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACTR; AIB-1; AIB1; CAGH16; CTG26; KAT13B; RAC3; SRC-3; TNRC14; TNRC16; TRAM-1; bHLHe42; pCIP; nuclear receptor coactivator 3; CBP-interacting protein; amplified in breast cancer 1 protein; class E basic helix-loop-helix protein 42; receptor-associated coa		

For research use only.

GENE DESCRIPTION	nuclear receptor coactivator 3
GENE ID	8202
ORIGIN SPECIES	Human
ACCESSIONS	BC119001
VECTOR	pDONR223
NCBI SYMBOL	NCOA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr