

ZNF75A-zinc finger protein 75a Gene

Reference: PTXBC119652

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 75A		
For research use only.			

GENE DESCRIPTION	zinc finger protein 75a
GENE ID	7627
ORIGIN SPECIES	Human
ACCESSIONS	BC119652
VECTOR	pDONR223
NCBI SYMBOL	ZNF75A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtattttccaggagaatggagttattggtccctcagaaggccctctacaatgatgtaagcaggaaaaatgagactgtcatctctagcattgtttgtctcccaaacctaaagtgatctctgtctagagcaaggggaagagccatgggtcaagtatccccggagttaag
gatagtgcggaaaatctctacagggttaagctcaaaaacgacactgaaaatcatcagcctgtgtctctttctgacttagaatacaagcatcagcaggcgtcatatcaaaaaagccaaagttaaagtctccagaaaacagcaggcaagaaaatcatttgatagcacagagtg
aaaatggcaccagaattttccagtgaagaaaagaaacattcaacctggaacaagagctgtcaaaactatggatcgtcacaagaagattgtgcaagagagaagccttttaaatgtcagggaatgtgggaaaacctcagagtttagctctgaccttattaaacccaagaattcacact
gaagagaaacccctataaatgtcaacagtggtgataagaggttttagatggagttcagatcttaataagcacttaacaacacacacaggaataaaacataataatgttcatgtgtgggaaaagcttcagtcatacaaatttacatacacacacaaagaactacacaggagaaaagccctca
catgtcatgaatgtgaaaaaattcagtcagaactccaccttattaaacacccggagaacccacacaggtgagcagccataactgttagcatatgcaggagaaacttcagcaggcggtaagcctcttagacaccagaactccacctgtga