

CLC-Charcot-Leyden crystal protein Gene

Reference: PTXBC119711

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GAL10; Gal-10; LGALS10; LGALS10A; LPPL_HUMAN; galectin-10; Charcot-Leyden crystal protein; eosinophil lysophospholipase; lectin, galactoside-binding, soluble, 10; lysolecithin acylhydrolase; Charcot-Leyden crystal galectin		

For research use only.

GENE DESCRIPTION	Charcot-Leyden crystal protein
GENE ID	1178
ORIGIN SPECIES	Human
ACCESSIONS	BC119711
VECTOR	pDONR223
NCBI SYMBOL	CLC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccctgctacccgtgccatacacagaggctgaccttctgtctactgtgtctactgtgacaatcaaggcgaccactgtctgtttcttgaatgaaccatactgcaggtggattccacactgagatgaaggagggaatcagacattgtcttcattccaagtgtgcttggctgctgtgtggtc
atgaacagccgtgagtagtggggcctggaagcagcaggtggaaatccaagaacatgcccttcaggatggccaagaatttgaactgagcatctcagtgctgccagataagtaaccaggttaatgggtcaafggccaatcctcttacaccttggaccatagaatcaagcctgaggctgtgaagatgt
gcaagtgtggagagatatctccctgaccaaatttaatgtcagctatttaagagataa