

ATP6V0E1-ATPase, H⁺ transporting, lysosomal 9kDa, V0 subunit e1 Gene

Reference: PTXBC119714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6H; ATP6V0E; M9.2; Vma21; Vma21p; V-type proton ATPase subunit e 1; ATPase, H ⁺ transporting, lysosomal 9kDa, V0 subunit e1; H(+)-transporting two-sector ATPase, subunit H; V-ATPase 9.2 kDa membrane accessory protein; V-ATPase H subunit; V-ATPase M9.2 s		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 9kDa, V0 subunit e1
GENE ID	8992
ORIGIN SPECIES	Human
ACCESSIONS	BC119714
VECTOR	pDONR223
NCBI SYMBOL	ATP6V0E1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgatcacggcctcactgtgccttcattgtgatgagcgtgttctgggcttcgtggcttcttggtgccttggttcatccctaagggtcctaaccggggagttatcattaccatgttggtgacgtgttcagttgctgctatctcttttggctgattgcaattctggcccaactcaaccctctctt
ggaccgcaattgaaaaatgaaaccatctgtatctgaagtatcattggccttga