

C10orf40-chromosome 10 open reading frame 40 Gene

Reference: PTXBC119749

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 10 open reading frame 40		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 40
GENE ID	283025
ORIGIN SPECIES	Human
ACCESSIONS	BC119749
VECTOR	pDONR223
NCBI SYMBOL	C10orf40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcataatacaacacactgtcatgtacaaaaagaagaatcaactcagcttagacgggtttatagttcatcagccagaggttggcttattcttttttaaggacgtagtctataactacatcacagaaaaacaacatgtagggtgcccttactcaagaaaggccttgacactggag
gtaagtttgaagaatatggaatgtcaacctgagggaattgggaagatgcaaatatggatctgtcaacctgagggaattgggaagatgcaaatagataaaaaacagaactatggtccaacagatgtgtatatattggtatgggtgattagatacaccacaaattcttgcgccatcctagtaa
gaggcatggcctaa