

DUSP21-dual specificity phosphatase 21 Gene

Reference: PTXBC119755

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LMWDSP21; dual specificity protein phosphatase 21; BJ-HCC-26 tumor antigen; LMW-DSP21; low molecular weight dual specificity phosphatase 21; dual specificity phosphatase 21		

For research use only.

GENE DESCRIPTION	dual specificity phosphatase 21
GENE ID	63904
ORIGIN SPECIES	Human
ACCESSIONS	BC119755
VECTOR	pDONR223
NCBI SYMBOL	DUSP21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagcatccgcgtctcttttcatcatcaggggtgccagcagccctccatctacagcttcccaataaccagaagcttgtttctcagcaatgggtgtgccgccaacgacaaactccttctgtccagcaatcgcacccgcatgtcaatgcctcgttgaagtgtcaacgtattct
tcgaggcattcagtaacataaagggtgcctgttaccgatgctgctgactcgcgtctctacgactttttgaccccatgtgatcttatccacaccatcgatataggcaggccgtacgctgctgcatgcatggctggagtgcggttccgctcactgtgccttgcgtacatgaaatacc
actccatgctcgtgctggagcccatatgacgaagtcgcgcgccccatccgccaacacggcttttgggaacagctcatcaattacgaattcaagctgtttaataacaacaccgtgcgcatgacaaactgccggtaggtaacatccctgacatctatgagaaggacctacgta
tgatgatataatgtaa