

LOC440356-hypothetical LOC440356 Gene

Reference: PTXBC119809

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDIPT antisense RNA 1 (head to head)		

For research use only.

GENE DESCRIPTION	hypothetical LOC440356
GENE ID	440356
ORIGIN SPECIES	Human
ACCESSIONS	BC119809
VECTOR	pDONR223
NCBI SYMBOL	LOC440356
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtccactccgatggggaatgaggggagagaagagcagctggccatctcaagctgcacctccttgagaggaggtccgcttcgttatctctgaggataacctgtccagatcagtcagaactcatggaggaggcttgtgactgcttgccttgaacctgtgcc
catcaaaaaaaaaagcagtcacaagaccaagcagctcagatatccaacgcgcattcttcaccaagacctga