

CST11-cystatin 11 Gene

Reference: PTXBC121080

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CST8L; CTES2; SC13; dJ322G13.6; cystatin-11; SC13delta; cystatin-related epididymal-specific protein; cystatin 11		

For research use only.

GENE DESCRIPTION	cystatin 11
GENE ID	140880
ORIGIN SPECIES	Human
ACCESSIONS	BC121080
VECTOR	pDONR223
NCBI SYMBOL	CST11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggctgagccctggcaggccctacagctcgtgtggccattctattgactctgatggccctccctaccaagcaaggagaaaaccttctaagcgtccatgaagtgtgagcagtagaaaactatgccaggacagcttcagtgatcaccgaccagtatacaaggaaagtgatga
caagtaccacttcaggatctccgagctctaaagtcagaggcaggtcactgaccacctggagtacacctgaacgtggaaatgcagtgaccacctgccaaaagccagagaccacgaactgtgtcccccaggaaaggagcttcacaagcaagtcactgcttctcagtggttgcctg
taccctggttgaacagtacaaaattttgaacaaaagctgcagcagtgactag