

KRTAP9-4-keratin associated protein 9-4 Gene

Reference: PTXBC121094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP9.4; KRTAP9.4; keratin-associated protein 9-4; keratin-associated protein 9.4; ultrahigh sulfur keratin-associated protein 9.4; keratin associated protein 9-4		

For research use only.

GENE DESCRIPTION	keratin associated protein 9-4
GENE ID	85280
ORIGIN SPECIES	Human
ACCESSIONS	BC121094
VECTOR	pDONR223
NCBI SYMBOL	KRTAP9-4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccactgtgtgtccctgtgtgtcagctacatgtgtcaggaccactgtgtcaggaccactgtgtgaagccaccactgtgtgaccactgtgcagcagcacacctgtgtccagccctctgtgtgtgtccagctgtgtccagctgtgtccgccaactgtgtcaaaacactg
ctgccagcccaactgtgtgaccagctgtgtccagcttctgtgtgcagcacacctgtgtgccagcccaactgtgtgtgtccagctgtgtgaccagcagctctgtgtgcacctgtgtactgcagaagaacctgtactacccacaactgtgtgcctgtgtgtgcttaaccagagctgtg
gtccaactgtgtccagccctgtgtgtccagccagctgtgtgagaccactgtgtccagcccaactgtgtgtacagctgtgtcagcctttttgtgtgtga