



Reference: PTXBC125083

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HD10; histone deacetylase 10		
For research use only.			

GENE DESCRIPTION	histone deacetylase 10
GENE ID	83933
ORIGIN SPECIES	Human
ACCESSIONS	BC125083
VECTOR	pDONR223
NCBI SYMBOL	HDAC10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr