



Reference: PTXBC125095

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HH17; protein sprouty homolog 4; sprouty homolog 4; sprouty RTK signaling antagonist 4		

For research use only.

GENE DESCRIPTION	sprouty homolog 4 (Drosophila)
GENE ID	81848
ORIGIN SPECIES	Human
ACCESSIONS	BC125095
VECTOR	pDONR223
NCBI SYMBOL	SPRY4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1