



Reference: PTXBC125113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf153; DJ964F7.1; disintegrin and metalloproteinase domain-containing protein 33; ADAM 33; a disintegrin and metalloprotease 33; a disintegrin and metalloproteinase domain 33; disintegrin and reprolysin metalloproteinase family protein; ADAM metallope		

For research use only.

GENE DESCRIPTION	ADAM metallopeptidase domain 33
GENE ID	80332
ORIGIN SPECIES	Human
ACCESSIONS	BC125113
VECTOR	pDONR223
NCBI SYMBOL	ADAM33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr