

GDEP-gene differentially expressed in prostate Gene

Reference: PTXBC125139

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GDEP; PCA4; PCAN1; prostate cancer associated transcript 4 (non-protein coding)		
For research use only.			

GENE DESCRIPTION	gene differentially expressed in prostate
GENE ID	118425
ORIGIN SPECIES	Human
ACCESSIONS	BC125139
VECTOR	pDONR223
NCBI SYMBOL	GDEP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgctgtgaaatctactacggtttgctgttttgaataaggagaaaaagagtgaggaaactgagaacatggatggccttgggaacgtggaaaagggtcactga