

CLDN23-claudin 23 Gene

Reference: PTXBC125148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLDNL; hCG1646163; claudin-23; 2310014B08Rik; claudin 23		
For research use only.			

GENE DESCRIPTION	claudin 23
GENE ID	137075
ORIGIN SPECIES	Human
ACCESSIONS	BC125148
VECTOR	pDONR223
NCBI SYMBOL	CLDN23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggacgccgggtgatgacgtgggcatggtgtggcgccctgcccctctgtcaacctgacggcaccctggcgccggctggcggtggtgaaggcttctgaaccagccagtgagctggagttgtaccagggcctgtggacatgtgtcgcgagcagacgacg
cgagcgcgagtgccgagacggaccagtggggctactcagggccagcccgtgctgtggcgccggcactcatgtcacctcgtggcgccacggctcctgggcttctgtggcgtcgtggcggtgcgctgtggcagacgagcccaactcgtgtggcagggctctcgg
ggctcgtgctctcgtcgtggcctcctcggcctcatcccggtgtcctgtgtacaaccttcttgggggacgcgacgtgtgctgcccggcccgccagccgggtcacgtgaggtcagctacagcctggtcctgggctacctgggcagctgcccctcgtcgtggcggttctcgtggtg
gtcagcttcgcccgtggcgagcagcgtgtgtcggcgccgcgcaaggacccctccggggcctcgcgcgacgagcgtcagcaccatcaaagtgagtgcccgagcccgacgtggcgccgccaatcaagtactacagcagcggccagcaccgacccgcccgtgcccagca
ccgcaagcccaagcccaagcccaaggtcggcttccccatgccgcccgcgcgcccgaaggcctacaccaactcgttgacgtcctcgtacggggaggggtgggagtgccagcgtccctcgtgcagcaccaccctgcgacagctcgtgcccctgcgactccgactctag