



Reference: PTXBC125176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	acyl-coenzyme A synthetase ACSM2A, mitochondrial; A-923A4.1; ACSM2; Homolog of rat kidney-specific (KS); acyl-CoA synthetase medium-chain family member 2; butyrate--CoA ligase 2A; butyryl-coenzyme A synthetase 2A; middle-chain acyl-CoA synthetase 2A; acyl		

For research use only.

GENE DESCRIPTION	acyl-CoA synthetase medium-chain family member 2A
GENE ID	123876
ORIGIN SPECIES	Human
ACCESSIONS	BC125176
VECTOR	pDONR223
NCBI SYMBOL	ACSM2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr