



Reference: PTXBC125269

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf84; NCS2; UPF0432; cytoplasmic tRNA 2-thiolation protein 2; cytosolic thiouridylase subunit 2 homolog; cytosolic thiouridylase subunit 2		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 84
GENE ID	348180
ORIGIN SPECIES	Human
ACCESSIONS	BC125269
VECTOR	pDONR223
NCBI SYMBOL	C16orf84
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]