



Reference: PTXBC126116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BGCAN; betaglycan; transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGF-beta receptor type III; TGFR-3; betaglycan proteoglycan; transforming growth factor beta receptor III; transforming growth factor beta receptor 3		

For research use only.

GENE DESCRIPTION	transforming growth factor, beta receptor III
GENE ID	7049
ORIGIN SPECIES	Human
ACCESSIONS	BC126116
VECTOR	pDONR223
NCBI SYMBOL	TGFBR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr