



Reference: PTXBC126117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTP-H1; PTPH1; tyrosine-protein phosphatase non-receptor type 3; cytoskeletal-associated protein tyrosine phosphatase; protein-tyrosine phosphatase H1; protein tyrosine phosphatase, non-receptor type 3		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 3
GENE ID	5774
ORIGIN SPECIES	Human
ACCESSIONS	BC126117
VECTOR	pDONR223
NCBI SYMBOL	PTPN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr