

C2orf73-chromosome 2 open reading frame 73 Gene

Reference: PTXBC126246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C2orf73; chromosome 2 open reading frame 73		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 73
GENE ID	129852
ORIGIN SPECIES	Human
ACCESSIONS	BC126246
VECTOR	pDONR223
NCBI SYMBOL	C2orf73
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaaagcctctgtggaatagtagcacttgccctcctgtgtacatcagcagaacttcaaaacaattttatagaatacatctctttatatacaatgatgccagaaaaactccgaatgagcctctccagggaagagacacggagcctttgtacagagagagataaaaccaggcagtagg
ccaacagttcctaaaggagcagagggtattactgaacactccagggtcacgttcacagacagtcacaaaaaacagagaaagaaactcagcggaaagcagaatgatctcaccagggtctctgccacaataatctcaggagctgttagagcctaaactcacttatcagaacagacgtca
gacaggcagccaaggcatgccccagcactcctgagagcagagaaaagacctcaggcgccactcaaaactgtaggagatgctctttcactaggcacaagcctttaaatccaccaattaaaaaatcagaataa