



Reference: PTXBC126354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PEK; WRS; eukaryotic translation initiation factor 2-alpha kinase 3; PRKR-like endoplasmic reticulum kinase; pancreatic EIF2-alpha kinase; eukaryotic translation initiation factor 2 alpha kinase 3		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 2-alpha kinase 3
GENE ID	9451
ORIGIN SPECIES	Human
ACCESSIONS	BC126354
VECTOR	pDONR223
NCBI SYMBOL	EIF2AK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr