

LNP1-leukemia NUP98 fusion partner 1 Gene

Reference: PTXBC126362

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NP3; leukemia NUP98 fusion partner 1		

For research use only.

GENE DESCRIPTION	leukemia NUP98 fusion partner 1
GENE ID	348801
ORIGIN SPECIES	Human
ACCESSIONS	BC126362
VECTOR	pDONR223
NCBI SYMBOL	LNP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcacaaagatgatgatgatgatgtgtcttttcccaatggatgagcagcttctggggccacagctggagagaggagatcagagaggactcgggaacgccagctgcaagccacagtcacaggaaaacctccctgccctgccacttctgtgcttccagaattcca
tcctctgactgcatcctagaaggcattctcatgaggaccaagaatttcgatccgtagccacgtacgggattacagaaaatactcagaggatgggtcattcaaggagccactggaatcaaaaggaagatccattccaaattgagaattttcagagtccttgaacggcaactgtgcttta
gaaccaagcgttctgctcttgggacctgaaagcagaaggagagaatgaaagagaatgctgaggtggagataaaatccgaaagaaagtagaggagaagaggagctctaggaaagaagagcagtgagaagcacacatggctccctgtttgaaaaagggcctgaataa