

LOC441208-hypothetical gene supported by AK094370 Gene

Reference: PTXBC126382

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical gene supported by AK094370
GENE ID	441208
ORIGIN SPECIES	Human
ACCESSIONS	BC126382
VECTOR	pDONR223
NCBI SYMBOL	LOC441208
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaccgaaaccagcggcagtcgacggccacctgagtcgccccctctgctggagccagcgggggtgcctgcagccgggacacctctctcgcgtctctctcccgcccgctcagccgccagccttgccaccgccccgagaagagcgcgcggcgccg
actgccccctggggcgccgagcggcgccctggacgtgccccggcctctctggccggccgcccgcctcgcctcctagctcccgcttgcctctaccctcttggctctcggggcacagcgcgccggtccggagcaggcgccgacatgagcgcagcagcgcg
gcctggcccgccgtaacggcgacgcgcgctactcggctcggatctaccttcagtag