

RLN2-relaxin 2 Gene

Reference: PTXBC126415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2-RLX; RLXH2; bA12D24.1.1; bA12D24.1.2; prorelaxin H2; H2-preprorelaxin; relaxin H2; relaxin, ovarian, of pregnancy; relaxin 2		

For research use only.

GENE DESCRIPTION	relaxin 2
GENE ID	6019
ORIGIN SPECIES	Human
ACCESSIONS	BC126415
VECTOR	pDONR223
NCBI SYMBOL	RLN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctcgccgtgtttttccacctgctaggagctgttttactgaaccaattttccagagcagtcgaggactcatggatggaggaagtattataattatcgggcgcgaattagttcgcgcgagattgccattgctggcatgagcacctggagcaaaaggctctgagccagggaagatgct
cctcagacacctagaccagtggcagaaattgtccatcctcatcaacaagatacagaaccataaatgatgtcagaattgttgctaattgccacaggagctgaagttaacctgtctgagatgcagccagcattaccacagctacaacaacatgtacgtgtattaaaagattccagctt
ctcttgaagaatttaagaaccttattcgcaatagacaaagtgaagccgcagacagcagtcctcagaattaaaatacttaggccttggatactcattctcgaaaaagagacaactctacagtgcatggcgaataaatgttgccatgttggtgtaccaaaagatctctgtagattttgctga