

INSL6-insulin-like 6 Gene

Reference: PTXBC126473

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	insulin-like peptide INSL6; RIF1; insulin-like peptide 5; insulin-like peptide 6; relaxin/insulin-like factor 1; insulin like 6		

For research use only.

GENE DESCRIPTION	insulin-like 6
GENE ID	11172
ORIGIN SPECIES	Human
ACCESSIONS	BC126473
VECTOR	pDONR223
NCBI SYMBOL	INSL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcggctcctcgttgcctgtgtggttgactcgtgtgttcgtgtgaactgagcgacatcagcagtgccaggaagctgtgcggcaggtacttggtaaagaaatagaaaaactctcggccatgccaaactggagccagttcgttcgaggaggaaaccccttc
tcacgggtgattgcacaggcctcggaaggtcgaaagctacagccatcaccagttcgaagcccgcaaacctctcccgccgggggaagaggcacaacccagtgctactcttgggaagcagtaaacagttgggaatgcagtcactacctgagtataaggataaaagg
gataatcacccttggttaagacaagagaattttctcatcacaataatcatgtatattcatgagaatgcaaaattcagaagaaacgtagaacacaaataaaacctaagcaattgttttggggcatcatcccaagaaaacgcagaggatattcagaaaagtgtgtcttacaggatga
caaaagaagaacttagcattgcattctccatatattgattttaaaggctaaaggaaaaagatcatcattgttaactaagataactaa