

TMEM107-transmembrane protein 107 Gene

Reference: PTXBC127649

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRVS638; PRO1268; transmembrane protein 107		

For research use only.

GENE DESCRIPTION	transmembrane protein 107
GENE ID	84314
ORIGIN SPECIES	Human
ACCESSIONS	BC127649
VECTOR	pDONR223
NCBI SYMBOL	TMEM107
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgggtctcagggcttggccctctcgcttctgacgtctcgtggcgcatctgggtgctcatcacttattctggtccgggacagcaacatacaggcctgcctgcctctcacgttcacccccgaggagtatgacaagcaggacattcagctgggtggccgcgtctctgtcaccctg
ggcctctttgcagtgagctggccggttctctcaggagctccatgttcaacagcaccagagcctcatctccattgggctcactgtagtgcacccgtggccctgtccttctcatattcgagcgttgggagtgactacgtattggatattttgtctctgcagtgcccttcagctgtcact
gaaatggctttattcgtcaccgtctttgggctgaaaaagaaccctctga