

MGC26718-similar to ankyrin repeat domain 20A Gene

Reference: PTXBC127726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 20A-like protein MGC26718; Putative ankyrin repeat domain-containing protein 20A-like protein MGC26718; MGC26718; ANKRD20A5; ankyrin repeat domain 20 family, member A5, pseudogene		

For research use only.

GENE DESCRIPTION	similar to ankyrin repeat domain 20A
GENE ID	440482
ORIGIN SPECIES	Human
ACCESSIONS	BC127726
VECTOR	pDONR223
NCBI SYMBOL	MGC26718
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttgttcggttcaggagccgcagggccagacggctcctgggtccatagaccacctctacacgggttcgggtaccgaatccggtactcggaaactgcagaagatccacaaggcagctgtcaagggcgacgtgcggagatggagcgtgtctggcgccaggagcggaga
ccttgacgccctggacaagcagcacagaactgctctacattggcctgtgccagtggccatgtgaaagtggctactctcctggttaacagaaaaatgacagttgatatctatgacaaagaaaatagaacgccttgatacaggtgtccattgccagggaagaggctgtgcccgttattctgctgg
aacatgggtccaatcaaaccttaaggatattacggcaacactgctctcattatgctgtgtatagtgagagcacctcactggcagaaaaactgctttccatggtgaaaatattgaagcactggacaaggtatag