

LOC90835-hypothetical protein LOC90835 Gene

Reference: PTXBC128129

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC90835
GENE ID	90835
ORIGIN SPECIES	Human
ACCESSIONS	BC128129
VECTOR	pDONR223
NCBI SYMBOL	LOC90835
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaaccgcaagaccagccattttgggaatgaggggtgcagtcggagcttgaacatctctccgagctgcggcggaagcgggggaagcgcagctcagtgcatgggtcggcgccgacgcgtgcgagtgctgggactcagtgagcagcggcgccggcggaagcg
gatgaagaccccgagccaactgtttccgcccgtgccccgacccggatctgcatgtggaagtacgtggagctccattccatgcaccagctggagaagaccaccaatctgagatgaggaggtgctggctgagctgctggagtaggtgtcctgagcagagcctgaggac
gcatcaccctggacctcttctgccacgcgtcattttctgccgacagggcttctcactggagcagacgtcagcggcttctgcccctgctccaggtcttcacaagctgtatgtcaaccccttgggcaacgtggagagtgctaccgctacttcaccagtgtcctttttgcatggagt
caggcgccctcttgcagctcagctcttcaagagggaacagctgtggccttgggaagactacgtgtgtaacacttacttccgcaactcaagctctataatagcttctcacaccccgagtgccgctggtatctgttgaacttactgggctacagccaccaactgtggcagagagt
gagacggagaagaagaagcaaggagatggaggagcaggcagttaccccgagaaagaggaaactagacagtggtcccccacagagccagagcccaatccacgtctccgagcctacatcaagccaagtgaacaaagagctggagcagctccagggcgtgtgg
aggagcggctcaaggccagcgaggaaggtcagcagcaagttgactgactagagcggccttccagctacttccggtaaggcaagagcaagaccaagtga