

## **TBC1D27-TBC1 domain family, member 27 Gene**

Reference: PTXBC128223

|                |                              |                     |                |
|----------------|------------------------------|---------------------|----------------|
| Product type:  | DNA                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | TBC1 domain family member 27 |                     |                |

**For research use only.**

---

|                         |   |
|-------------------------|---|
| <b>GENE DESCRIPTION</b> | TBC1 domain family, member 27   |
| <b>GENE ID</b>          | 96597   |
| <b>ORIGIN SPECIES</b>   | Human   |
| <b>ACCESSIONS</b>       | BC128223  |
| <b>VECTOR</b>           | pDONR223  |
| <b>NCBI SYMBOL</b>      | TBC1D27   |
| <b>LENGTH</b>           | ProteoGenix bp  |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### **ORF (Only for reference, product provided without STOP codon )**

atggagatggatgagcaccggatctctgccaccaggggcaaggcaacatcatcattactaagtatgagcagggacgctga