

KIAA0355-KIAA0355 Gene

Reference: PTXBC128232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA0355		
For research use only.			

GENE DESCRIPTION	KIAA0355
GENE ID	9710
ORIGIN SPECIES	Human
ACCESSIONS	BC128232
VECTOR	pDONR223
NCBI SYMBOL	KIAA0355
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgagcataatctcctgaatccaaaggcgattgaagcaagtttgcagggctgtgcagcagggcggaagcccgagcagcagggcgaggcagacaccccgagcccatgcagtgtagctccccaccgtccctgtgcagataggatcgacttctgaaggcgctctctta
atgagtcggccgacaaatctgaaactaagacgcatacaatgttacagctgatgaaggaggcagggctgtataatggaatcacatccaggatgatttctctgtgactgaagtgtgaacaggtttgccctccacatggcgaggtgcctgcaagacggcggtgcagctgtgttggc
caggctggactggtgtgtgtgacacagcacagattgagaataaagaagcctatgccccagatcagtttagaaggctctagaatcgtggtcaagtcctccatccatgggtgcctgaagaagacccctgctaccatgtccctgctgcagagaagccttgatcctgagaagaccctgggtct
agtgagcgtgctctacacagctgtgctggacctaaccgctggaggctggaaggagcaagctctctctctttaa