



Reference: PTXBC128260

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	vitrin		
For research use only.			

GENE DESCRIPTION	vitrin
GENE ID	5212
ORIGIN SPECIES	Human
ACCESSIONS	BC128260
VECTOR	pDONR223
NCBI SYMBOL	VIT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr