

KRTAP23-1-keratin associated protein 23-1 Gene

Reference: PTXBC128435

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP23.1; keratin-associated protein 23-1; keratin associated protein 23-1		
For research use only.			

GENE DESCRIPTION	keratin associated protein 23-1
GENE ID	337963
ORIGIN SPECIES	Human
ACCESSIONS	BC128435
VECTOR	pDONR223
NCBI SYMBOL	KRTAP23-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctacaactgctgctgtgaaacttctctccatctgtgagggctacctgtgtactcaggctactccgtgggtgctctctgtacccagcaacctggtctacagcactgaaccttggatctccagcacctgccagctgggtctctctctgcaagggcttccaggagactgtctgg
gaaaccctag