

DEFB105A-defensin, beta 105A Gene

Reference: PTXBC128437

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD-5; DEFB-5; DEFB105; beta-defensin 105; beta defensin 5; defensin, beta 5; defensin beta 105A		

For research use only.

GENE DESCRIPTION	defensin, beta 105A
GENE ID	245908
ORIGIN SPECIES	Human
ACCESSIONS	BC128437
VECTOR	pDONR223
NCBI SYMBOL	DEFB105A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgatcagggaagacattttatttctatttctatgtttcttcaatttgggtcaactgccatcagggtgccaggcaggacttgattttcccaaccatttcacatcagggtgagtttgcgtctgtgagtcgtcaagcttggtcgggaaaaatgcagggaaggagtgcttgagaatgagaagccc
gatggaaatgcaggctgaactttctctgctgcagacagagatctga