



Reference: PTXBC129807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-ALPHA3; protocadherin alpha-3; KIAA0345-like 11; PCDH-alpha-3; protocadherin alpha 3		

For research use only.

GENE DESCRIPTION	protocadherin alpha 3
GENE ID	56145
ORIGIN SPECIES	Human
ACCESSIONS	BC129807
VECTOR	pDONR223
NCBI SYMBOL	PCDHA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]