



Reference: PTXBC130283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	peptidyl-cysteine S-nitrosylase NOS2; HEP-NOS; NOS; NOS2A; nitric oxide synthase, inducible; NOS, type II; hepatocyte NOS; inducible NO synthase; inducible NOS; nitric oxide synthase 2, inducible; nitric oxide synthase 2A (inducible, hepatocytes); nitric		

For research use only.

GENE DESCRIPTION	nitric oxide synthase 2, inducible
GENE ID	4843
ORIGIN SPECIES	Human
ACCESSIONS	BC130283
VECTOR	pDONR223
NCBI SYMBOL	NOS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr