



Reference: PTXBC130293

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P76; transmembrane 9 superfamily member 2; 76 kDa membrane protein; dinucleotide oxidase disulfide thiol exchanger 3 superfamily member 2; transmembrane protein 9 superfamily member 2		

For research use only.

GENE DESCRIPTION	transmembrane 9 superfamily member 2
GENE ID	9375
ORIGIN SPECIES	Human
ACCESSIONS	BC130293
VECTOR	pDONR223
NCBI SYMBOL	TM9SF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]