

LOC440337-hypothetical LOC440337 Gene

Reference: PTXBC130355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC440337
GENE ID	440337
ORIGIN SPECIES	Human
ACCESSIONS	BC130355
VECTOR	pDONR223
NCBI SYMBOL	LOC440337
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctgtctctgacaaatcctggtgcttcacgtcctttaaatacaaacgaggtgatctagtgttctctctgctgctctgtgcttctagcatcctggatagccgtgtatgtttagaagaggggtccaccctgtgccccctaaagatgatggctcattatcaagaccaggctggag
gctcagctaccagagcttctgttttggagggaaggctgcagcagcagatgttgcctcgcctcagctccctgcagccttggagctccatctgcccacactcgggcacagctgccagtgcttccccctccgtgctcacccttccactcctcagcagagctgggtttcatgcatctgggccc
aacagagagatggtgaccagaccctgcgacgtggtgccacaatcctccaacagctgctgggcatcttctgcagtagcctgtcagcatgctggaagcagcaccgatggaggtgttga