



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

VN1R2-vomeronasal 1 receptor 2 Gene

Reference: PTXBC130356

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

V1RL2; vomeronasal type-1 receptor 2; G-protein coupled receptor GPCR25; V1R-like 2; V1r-like receptor 2; hGPCR25; pheromone receptor; vomeronasal 1 receptor 2

For research use only.

GENE DESCRIPTION	vomeronasal 1 receptor 2
GENE ID	317701
ORIGIN SPECIES	Human
ACCESSIONS	BC130356
VECTOR	pDONR223
NCBI SYMBOL	VN1R2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr