

PAGE2B-P antigen family, member 2B Gene

Reference: PTXBC130359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT16.5; GAGEE3; PAGE-2B; G antigen, family E, 3; P antigen family, member 2B; prostate-associated gene 2B protein; PAGE family member 2B		

For research use only.

GENE DESCRIPTION	P antigen family, member 2B
GENE ID	389860
ORIGIN SPECIES	Human
ACCESSIONS	BC130359
VECTOR	pDONR223
NCBI SYMBOL	PAGE2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagcatgtgagaacaagatccaatcctcagaaagaggaaatgaccaagagctctccagccagttggatctgtgattgtccaggagccactgaggaaaaacgtcaagaagaggaaaccacactgataatcagggtattgcacctagtggggagatgaaaatgaaggagc
acctgccgttcaaggcctgacatggaagctttcaacaggaaactgctctgcttaagatagagatgagcctggagatggtcctgatgtcaggaggaggattatgccactttgatctactaaagtgcctggaagcagggtgatgcgaaccatag