



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C10orf25-chromosome 10 open reading frame 25 Gene

Reference: PTXBC130369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf25; chromosome 10 open reading frame 25		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 25
GENE ID	220979
ORIGIN SPECIES	Human
ACCESSIONS	BC130369
VECTOR	pDONR223
NCBI SYMBOL	C10orf25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccggggccaccgaaagcggtgtgttcttctgtgttcttcttctcccccaactcggaaagcgctgtgtgaccccgggagactcaagagctgtgaaccgccctgtgtgtgtcaggcgtttgaaacaaactctagctctcaacctagaactcgtacattccacaaa
atttgagattttgagaccctgtgactctttcacgaagtcacctaattatcacaaaagaagctctgtcaagatggagatgtgtatactcaattttagagaaaggaacaaagcctcagagaagaagtgctctcgttaccgtgtgtcctctcctcagccactctgccgtgtcctcagtcactga