



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PRB4-proline-rich protein BstNI subfamily 4 Gene

Reference: PTXBC130386

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	basic salivary proline-rich protein 4; BstNI type basic salivary proline-rich protein 4; parotid o protein; salivary proline-rich protein II-1; salivary proline-rich protein Po; proline rich protein BstNI subfamily 4		

For research use only.

GENE DESCRIPTION	proline-rich protein BstNI subfamily 4
GENE ID	5545
ORIGIN SPECIES	Human
ACCESSIONS	BC130386
VECTOR	pDONR223
NCBI SYMBOL	PRB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgctgattctctgctcagtgccctctgctgcccctgagctcagctcagagatgctcaagtgaaagatcagccaggaagaatctctctctaatafatcagaagccagaaggacgacgccacaaggaagaaacacgcccacctctccaggaaagccacaaggacc
acccccacaaggaggaaacagatcccaaggctcccccactctctccaggaaagccagaaggacgacccccacaaggaggcacaacgattccaaaggtccccacctctatccaggaaagccagaagaagcacccccacaaggaggaaacacagttcccaaggtaccccactctctcag
gaaagccagaagacgacccccacaaggaggcaacaacagttccaaaggtccccacctctatccaggaaagccagaagaagcacccccacaaggaggcaaacaggtccccacctctctcaggaaagccagaagaagcacccccacaaggagggaacacagttcccaaggt
tccccacctctatccaggaaagccagaaggacccccacagaggaagaaacacagttccaaaggtccccacctctctcaggaaagccagaagacccccacaagaagaagcaacaagctctcaaggtccccacctctgaaagccacaaggcccacccaccagaggga
atccccagcctcagctcagccactctcgtgaaagcccccagggccacccactctcccaaggggctcagcgcacccacagacgtcccaaggagacaacacgctctcccaagtaa

More info: Please contact us for any question: contact@proteogenix.fr