



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

MPN2-marapsin 2 Gene

Reference: PTXBC130400

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine protease MPN2; MPN2; serine protease 38; marapsin 2; protease, serine 38		
For research use only.			

GENE DESCRIPTION	marapsin 2
GENE ID	339501
ORIGIN SPECIES	Human
ACCESSIONS	BC130400
VECTOR	pDONR223
NCBI SYMBOL	MPN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcctccctgtcttcgcatgggcccactcgggcccctctgccggcgctctctctgctgctgctctgggtggcgccctcccccgggtgcgacatgtgtccacagacagccagagaaacagggaatctccctaactggcagcgtggcgctgtgtgctgcccagcatggaggggaaat
ctggggcggcgtctctcgtgccccagaggaagtggcgctgctgcagctcagctgcacacacaggcctccacgtctccatctcaatgaatctatggctgtgctgctgcagctgctgcagcgcctgtctccagggacaagaatacaaaattatgatactctatgctagctgcctgtaaact
cagggtgtggacggcaaacacacacagtggtatgaggtgaacagggtgatctctgcacccacatagatgatattacaccccactcggagtgacgtgctggcctgtgcagctgaagaccgcatctgtttttctgactgcgtccctgttttgcctcaactcagaagtgaaacttaacct
gcaattgtcagctcgggatgggcatgctcctgaatgctcaaaacaaggctgacacctcagacgagctcaggaggtgcagctccctgctgacttcaggagccctggcgtgcacactgctctacgacacattctctacatcatgccgacatctgtgtgctgtgggacatctgaatctgaagaccgtg
tctggggcgactccggggcccactctgtctgtaactcaaccgcagctgctgctgcagattgaattgtgagctgcggccgagcgtctctcaccccctctgacactggaggtgatgcaggttttctatctcaaaatggatattgataacatagaatacacgccactctctcagccagc
ccctctctctctcaactctctggcccactcgaactctctcaatctgacccactctgctctgctcactctgaa