



Reference: PTXBC130409

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAGI-3; dJ730K3.2; membrane-associated guanylate kinase, WW and PDZ domain-containing protein 3; membrane-associated guanylate kinase inverted 3; membrane-associated guanylate kinase-related 3; membrane associated guanylate kinase, WW and PDZ domain conta		

For research use only.

GENE DESCRIPTION	membrane associated guanylate kinase, WW and PDZ domain containing 3
GENE ID	260425
ORIGIN SPECIES	Human
ACCESSIONS	BC130409
VECTOR	pDONR223
NCBI SYMBOL	MAGI3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr