

LOC284297-hypothetical LOC284297 Gene

Reference: PTXBC130538

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S5D-SRCRB; soluble scavenger receptor cysteine-rich domain-containing protein SSC5D; scavenger receptor cysteine rich domain containing (5 domains); scavenger receptor cysteine rich family, 5 domains; scavenger receptor cysteine-rich glycoprotein; soluble		

For research use only.

GENE DESCRIPTION	hypothetical LOC284297
GENE ID	284297
ORIGIN SPECIES	Human
ACCESSIONS	BC130538
VECTOR	pDONR223
NCBI SYMBOL	LOC284297
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgagccacctgccctgtgagctggtgctgtgtgagggatgtgggtggtcagctgcagagactgacccaggtctgtggaacaggagcggcaggagcggccaagccctgctgtgggctgacgcagctggtagaagctgccggggtctgggcagctgggtgagg
ctgtgaagagactggcagagatggcctggaccaccagcatgcctgcaccaaccaccactacccagaggagaagaagaagacccctgaggggagacgtgtga