

FLJ26850-FLJ26850 protein Gene

Reference: PTXBC130542

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLJ26850 protein		
For research use only.			

GENE DESCRIPTION	FLJ26850 protein
GENE ID	400710
ORIGIN SPECIES	Human
ACCESSIONS	BC130542
VECTOR	pDONR223
NCBI SYMBOL	FLJ26850
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcactagacagtggggccctcaaggacctggcaaatggggcttggacctggatcggtccacgcgaagcttgcctccaccacgacatccccacgttggfggcgacgtgccccggccccacggatactccgcgcctgtcagactccctgatgaactaccctccagagtagc
cgccggagctcgggctctgagggcgacggctcctgatggcagatgcgggagaaactctggcgtcaggcggccctcgcgtggagcacacgaagctgtggttattctgcttcagtagtggggfggagaaggcgatccacgcagctgcgtctattctgtggatcaatcgaaaaata
cgttctgtaa