

C7orf59-chromosome 7 open reading frame 59 Gene

Reference: PTXBC130553

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0539 protein C7orf59; C7orf59; regulator complex protein LAMTOR4; late endosomal/lysosomal adaptor and MAPK and MTOR activator 4; late endosomal/lysosomal adaptor, MAPK and MTOR activator 4		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 59
GENE ID	389541
ORIGIN SPECIES	Human
ACCESSIONS	BC130553
VECTOR	pDONR223
NCBI SYMBOL	C7orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacttctgcaactgacccaggggctggagcgaatcccagaccagctcggctacctggtactgagtgaaggtgcagtgctggcgtcatctggggacctggagaatgatgagcaggcagccagtgcctctctgagctggtcagcacagcctcggtttccggtgcaccgcccgtgaa
tftgcccttcaagcgccctgtctgtggtctttggagaacacacactgctggtgacggtgtcaggacagaggggtgtttgtggtgaagaggcagaaccgaggtcgggagccattgatgtctga